

P/N 980-04926 Article Structure Description	
Vial type,	1.5ml Screw Neck Vial ND9 Kit includes Vials, Caps and Septa
Size	32 x 11.6mm
Glass material	Clear glass
Glass quality	1st hydrol. class
Volume	1.5 ml
Cap	Blue PP Short Thread Cap
Cap center hole size	6mm
Septum type	Silicone white/PTFE red ultra clean
Septum thickness	1.0 mm
Septum hardness	55° shore A
Temperature resistant	-60 to 200°C
Additional remarks	Vial with label and filling lines
Packing	100 vials with caps and septa in PP box

