

P/N 980-04988 Article Structure Description	
Vial type,	4.0 ml Screw Neck Vial ND13 Kit includes Vials, Caps and Septa
Size	45 x 14.7mm
Glass material	Clear glass
Glass quality	1st hydrol. class
Volume	4.0 ml
Cap	Black Polypropylene Screw Cap
Cap center hole size	8.5mm
Septum type	Silicone dark blue/PTFE white
Septum thickness	1.3mm
Septum hardness	45° shore A
Temperature resistant	-60 to 200°C
Additional remarks	
Packing	100 vials with caps and septa in PP box

