

P/N 980-08338 Article Structure Description	
Cap	Blue PP Short Thread Cap
Cap center hole size	6mm
Septum type	Silicone white/PTFE red ultra clean
Septum thickness	1.0 mm
Septum hardness	55° shore A
Temperature resistant	-60 to 200°C
Additional remarks	Cap fitting for Vial 980-08337. Same cap and septum as included in kit P/N 980-04926. The UltraClean Septum is recommended for trace, complex, or critical analysis.
Packing	1000 caps and septa in PP box

